

Store at
4°C

Sample Diluent A



Cell Signaling
TECHNOLOGY®

#71637

125 mL

Support: +1-978-867-2388 (U.S.)
www.cellsignal.com/support

Orders: 877-616-2355 (U.S.)
orders@cellsignal.com

New 03/21

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Sample Diluent A is specifically formulated to dilute serum and plasma samples for use with several of the Cell Signaling Technology (CST) SARS-CoV-2 serological and blocking antibody detection ELISA kits (#20154, #37322, #58873, #88005, and #69486). Serum or plasma samples should be diluted with Sample Diluent A as indicated in each kit's specific protocol.

Storage: Store at 4°C.

For product specific protocols and a complete listing of recommended companion products please see the product web page at www.cellsignal.com.

Thank you for your recent purchase. If you would like to provide a review visit www.cellsignal.com/comments.

www.cellsignal.com

© 2021 Cell Signaling Technology, Inc.

Cell Signaling Technology is a trademark of Cell Signaling Technology, Inc.

Applications: W—Western IP—Immunoprecipitation IHC—Immunohistochemistry ChIP—Chromatin Immunoprecipitation IF—Immunofluorescence F—Flow cytometry E-P—ELISA-Peptide **Species Cross-Reactivity:** H—human M—mouse R—rat Hm—hamster Mk—monkey Mi—mink C—chicken Dm—D. melanogaster X—Xenopus Z—zebrafish B—bovine Dg—dog Pg—pig Sc—S. cerevisiae Ce—C. elegans Hr—Horse **All**—all species expected Species enclosed in parentheses are predicted to react based on 100% homology.