

Store at Room
Temperature

#9808

Phosphate Buffered Saline (PBS-20X)



Supplied as a 20X stock

1000 ml

Orders ■ 877-616-CELL (2355)
orders@cellsignal.com

Support ■ 877-678-TECH (8324)
info@cellsignal.com

Web ■ www.cellsignal.com

rev. 11/07/17

For Research Use Only. Not For Use In Diagnostic Procedures.

Description: Phosphate Buffered Saline (PBS) formulated to be used as a buffer for ELISA. 1 X PBS contains 3.2 mM Na_2HPO_4 , 0.5 mM KH_2PO_4 , 1.3 mM KCl, 135 mM NaCl, pH 7.4.

Directions for Use: Dilute Phosphate Buffered Saline (PBS - 20X) to a 1X solution using ddH₂O. This product supplies enough 20X material to make 20 liters of 1X solution.

Storage: Supplied as a 20X solution. Store at room temperature.

Please visit www.cellsignal.com for a complete listing of recommended companion products.