

Store at Room Temperature  
**#9872**

# Phosphate Buffered Saline (PBS-1X) pH7.2 (Sterile)

- Small 10 ml  
 Large 50 ml



**Orders** ■ 877-616-CELL (2355)  
orders@cellsignal.com  
**Support** ■ 877-678-TECH (8324)  
info@cellsignal.com  
**Web** ■ www.cellsignal.com

rev. 11/07/17

**For Research Use Only. Not For Use In Diagnostic Procedures.**

**Description:** Phosphate Buffered Saline (PBS-1X) pH 7.2 (Sterile) is specifically formulated to be used as a buffer for reconstituting sterile cytokines or antibodies. 1X PBS contains 10 mM Na<sub>2</sub>HPO<sub>4</sub>, 10 mM NaH<sub>2</sub>PO<sub>4</sub>, 150 mM NaCl.

**Endotoxin:** Less than 0.1 EU/ml PBS.

**Storage:** Supplied as a 1X solution. Store at room temperature. Maintain sterility.